

How To Measure for the Looking Forward Head Support ®

The head 'LIFT' benefit is achieved by the length of the band and how it fits the recipient. For this reason, the band is made in different lengths so that there is enough 'Pull' but not too much "Play" to get the required head lift.

When the band or cap is correctly placed on the head (cap) or forehead (band), cross it behind the neck and then PULL THE ENDS OF BAND DOWNWARD TOWARD YOUR HIPS TO GET AS MUCH STRETCH FROM THE BAND AS POSSIBLE.

Keeping the band taut, bring the two ends of the band around to the front of your body and secure it under the breast. This pulling and stretching will assure you get maximum head lift. The crossing, pulling and stretching of the band is how the head lift is obtained.

The band should be worn on top of an undergarment (Tee Shirt or Camisole) so measuring should be done including the garment being worn with the Looking Forward.

Measurements are determined by measuring the patient around the torso, under the arms, just under the breast area.

Standard, in stock sizes are Extra Small - 24" -29" Small 30"- 36" Medium/Large 37"-42" Extra Large 48" +
Custom sizes are available by special order.

If the patient is measuring close to the end of the size selected, go up one size.

If the patient has postural stoop, allow enough length to accommodate the bent posture. If possible, it is helpful (but not mandatory) to measure the patient in a standing position for better accuracy.